

CROC®



Single Use Biopsy Forceps
For Gastroenterology Procedures

CHOOSE
CROC®
FOR THE
PERFECT
BITE!



primed

CROC® ... FOR THE PERFECT BITE!

CROC® Single Use Biopsy Forceps are engineered and manufactured to collect tissue samples of the upper digestive tract and colon for histological examination.

The devices are supplied sterile and are intended for a single use only.

Consider the following features of the CROC® and discover why it is the finest biopsy forceps available:

- Large capacity oval fenestrated jaws assure excellent quality tissue sampling for accurate histological diagnosis.
- Precise alignment of the fenestrated cups achieves a tight and even closure of the jaws.
- Ergonomic handle ensures rapid and controlled response to the actuation of the jaws.
- Resin coated sheath promotes a smooth insertion and passage of the forceps through the channel of the endoscope.
- Patented design and construction guarantees a frictionless actuation of the device even in the most tortuous anatomy.
- Jaws with or without a spike give the medical professional a preferred choice of forceps.

MODEL NUMBER	JAWS DIAMETER	LENGTH	JAWS TYPE
Gastroscope			
CR24G	2.4 mm	160 cm	No Spike
CR24GS	2.4 mm	160 cm	Spike
Colonoscope			
CR24C	2.4 mm	230 cm	No Spike
CR24CS	2.4 mm	230 cm	Spike

* Minimum channel size 2.8 mm

Prior to use read **Instructions for Use** pamphlet.

This product is protected by one or more of the following Patents:

U.S. Patents: 6,309,404; 7,171,839

Canadian Patent: 2,094,529

Other International Patents and Patents pending



Manufactured by:
PriMed Instruments Inc.
1080 Tristar Drive, Unit #14
Mississauga, ON L5T 1P1
Canada

Phone: 905-565-0565

Fax: 905-565-0566