

# PriMed Instruments Inc.

## ORDER FORM

### INSTRUCTIONS

#### STEP 1

Complete the **Shipping and Contact Information** section.

#### STEP 2

Indicate the quantity of product(s) you wish to order.

#### STEP 3

Save the completed order PDF to your computer.

#### STEP 4

Open the completed PDF and click on the **SUBMIT** button to open a new email page already addressed to PriMed's Ordering Center.

#### STEP 5

Attach the **ORDER FORM** PDF in the email and send.



## SHIPPING AND CONTACT INFORMATION

Hospital

Address

City

Province / State

Country

Postal / Zip Code

Telephone

FAX

Email

Your Name

Comments or Questions

Purchase Order Number

Delivery Date

## QUICK LINKS

To jump to a specific product group, click on one of the buttons below.

[Single Use CROC® Series Biopsy Forceps](#)

[Reusable BX Series Biopsy Forceps](#)

[Single Use Endoscope Cleaning Brushes](#)

[Reusable Endoscope Cleaning Brushes](#)

## Standard Capacity CROC® Single Use Biopsy Forceps

### For Gastroscopes

PRODUCT NUMBER	JAWS DIAMETER	WORKING LENGTH	OVAL FENESTRATED JAWS	MINIMUM CHANNEL SIZE	QUANTITY
CR24G-20	2.4mm	160cm	No Spike	2.8mm	
CR24GS-20	2.4mm	160cm	Spike	2.8mm	

### For Colonoscopes

CR24C-20	2.4mm	230cm	No Spike	2.8mm	
CR24CS-20	2.4mm	230cm	Spike	2.8mm	

*Note: Twenty devices in a dispenser box*

## High Capacity CROC® Single Use Biopsy Forceps

### For Gastroscopes

PRODUCT NUMBER	JAWS DIAMETER	WORKING LENGTH	ALLIGATOR DOUBLE FENESTRATED JAWS	MINIMUM CHANNEL SIZE	QUANTITY
H624S-20	2.4mm	160cm	Spike	2.8mm	
H624N-20	2.4mm	160cm	No Spike	2.8mm	
H628S-20	2.8mm	160cm	Spike	3.2mm	
H628N-20	2.8mm	160cm	No Spike	3.2mm	

### For Colonoscopes

H924S-20	2.4mm	230cm	Spike	2.8mm	
H924N-20	2.4mm	230cm	No Spike	2.8mm	
H928S-20	2.8mm	230cm	Spike	3.2mm	
H928N-20	2.8mm	230cm	No Spike	3.2mm	

*Note: Twenty devices in a dispenser box*

## Reusable BX Series Biopsy Forceps

### For Colonoscopes – Working Length 240cm

PRODUCT NUMBER	JAWS TYPE	JAWS DIAMETER	NEEDLE	QUANTITY
BX 100	Oval fenestrated	2.4mm	Yes	
BX 100N	Oval fenestrated	2.4mm		
BX 110	Alligator fenestrated	2.4mm	Yes	
BX 110N	Alligator fenestrated	2.4mm		
BX 140	Oval fenestrated	3.0mm	Yes	
BX 140N	Oval fenestrated	3.0mm		

### For Gastroscopes – Working Length 155cm

PRODUCT NUMBER	JAWS TYPE	JAWS DIAMETER	NEEDLE	QUANTITY
BX 200	Oval fenestrated	2.4mm	Yes	
BX 200N	Oval fenestrated	2.4mm		
BX 210	Alligator fenestrated	2.4mm	Yes	
BX 210N	Alligator fenestrated	2.4mm		
BX 300	Round fenestrated	2.4mm	Yes	
BX 300N	Round fenestrated	2.4mm		
BX 240	Oval fenestrated	3.0mm	Yes	
BX 240N	Oval fenestrated	3.0mm		
BX 420	Oval fenestrated	1.8mm		
BX 520	Alligator	1.8mm		

### For Sigmoidoscopes – Working Length 120cm

PRODUCT NUMBER	JAWS TYPE	JAWS DIAMETER	NEEDLE	QUANTITY
BX 250	Oval fenestrated	2.4mm	Yes	
BX 250N	Oval fenestrated	2.4mm		
BX 260	Alligator fenestrated	2.4mm	Yes	
BX 260N	Alligator fenestrated	2.4mm		
BX 280	Oval fenestrated	3.0mm	Yes	

**For Duodenoscopes – Working Length 180cm**

PRODUCT NUMBER	JAWS TYPE	JAWS DIAMETER	QUANTITY
BX 430	Oval fenestrated	1.8mm	

**For Bronchoscopes – Working Length 105cm**

PRODUCT NUMBER	JAWS TYPE	JAWS DIAMETER	QUANTITY
BX 410	Oval fenestrated	1.8mm	
BX 510	Alligator	1.8mm	

## Single Use Endoscope Cleaning Brushes

### Single-Use Channel Cleaning Brushes

PRODUCT NUMBER	BRUSH HEAD DIAMETER	OVERALL LENGTH	MIN. CHANNEL SIZE	COLOR	QUANTITY
9190	5-6.4mm taper	240cm	2.6mm	Navy Blue	
9192 Triple Head	3-6.3mm <b>DART®</b> Push and Pull	240cm	2mm	Navy Blue	
9115	3-5mm taper	115cm	2mm	Sky Blue	
<b>9025 MICRO BRUSH</b>	<b>2.5mm</b>	<b>120cm</b>	<b>1.2mm</b>	<b>Yellow</b>	<input type="text"/>
<b>9020 MICRO BRUSH</b>	<b>2mm</b>	<b>120cm</b>	<b>1.2mm</b>	<b>Orange</b>	<input type="text"/>
<b>9010 MICRO BRUSH</b>	<b>2.5mm</b>	<b>110cm</b>	<b>1mm</b>	<b>Twisted Wire</b>	<input type="text"/>

### Single-Use Suction Well and Valve Cleaning Brushes

PRODUCT NUMBER	BRUSH HEAD DIAMETER	WORKING LENGTH	HANDLE COLOR	QUANTITY
9990 Suction Well	11mm	5.5cm	Navy Blue	
9992 Double Ended Combo	4.3mm & 11mm	5.5cm & 5.5cm	Navy Blue	

### Single-Use Endoscope Cleaning Brushes

PRODUCT NUMBER	DESCRIPTION	QUANTITY
9199	Channel & Suction Well Cleaning Brush Kit (9190 + 9990) – 2 pieces	
9155	Channel & Suction Well/Valve Cleaning Brush Kit (9115 + 9992) – 2 pieces	

*Other sizes, custom configurations and combinations are available.  
Please contact us with your requirements.*

## Reusable Endoscope Cleaning Brushes

### Channel Cleaning Brushes

PRODUCT NUMBER	BRUSH HEAD DIAMETER	CATHETER DIAMETER	WORKING LENGTH	FOR USE IN:	QUANTITY
8107	3-5mm	1.6mm	70cm	Cystoscope	
8115	3-5mm	1.6mm	115cm	Bronchoscope	
8160	3-5mm	1.6mm	160cm	Slim Gastroscope/ Ped	
8170	6.5-7mm	1.8mm	240cm	Colonoscope/ Gastroscope	
8180	3-6mm	1.8mm	240cm	Colonoscope/ Gastroscope	
8182	Double End	1.8mm	240cm	Colonoscope/ Gastroscope	
8190	3-6mm	1.8mm	120cm	Sigmoidoscope	
8195	3-6mm	1.8mm	250cm	Colonoscope	

### Suction Well Cleaning Brush

PRODUCT NUMBER	BRUSH HEAD DIAMETER	WORKING LENGTH	QUANTITY
8880	11mm	5.5cm	

**SUBMIT**