

EZPZ-II Polyp Trap™

Easy To Use
2-Basket Polyp Trap

- ✓ Optically Clear Material
- ✓ Secure Polyp Capture
- ✓ Cost Effective



primedENDO.com



EZPZ-II Polyp Trap™

EZPZ-II Polyp Traps™ are engineered and manufactured to retrieve polyp tissue through the endoscope channel for histopathological diagnosis.

Devices are intended for single use only.



Consider the following features & benefits of the EZPZ-II Polyp Traps™

- ✓ **Optically clear material.**
 - Easy to see through to confirm presence of polyp tissue in baskets.
- ✓ **Two removable, selectable and numbered baskets.**
 - Economical evacuation solution for the most common number of polyps removed during colonoscopy.
 - Simple turning of the lid or the barrel selects the basket to be used.
 - Large numbers printed on baskets allow accurate identification of tissue.
 - Basket tabs provide for easy removal from the barrel.
- ✓ **Centrally located vacuum lumen.**
 - Fast connection of the device.
 - Baskets do not get knocked out during removal of the lid.
- ✓ **Fine diamond patterned strainer basket bottom.**
 - Even the smallest tissue samples are collected.
- ✓ **Optimal size barrel.**
 - Excellent flow of fluid without blockages.

PART NUMBER	DESCRIPTION	UNITS PER BOX
PLT001	Two Removable Baskets Polyp Trap	100

Read Instructions for Use pamphlet before using.



primedENDO.com

Manufactured by:
PriMed Instruments Inc.
1080 Tristar Drive, Unit #14
Mississauga, ON L5T 1P1
Canada

☎ 905-565-0565

☎ 905-565-0566

TO ORDER CONTACT

✉ sales@primedENDO.com

☎ Toll Free: 877-565-0565

☎ Toll Free: 877-565-0566

