



# MulticROC®

Multi-Sampling Biopsy Advantage

- Reduced Procedure Time
- Improved Diagnostics
- Excellent Tissue Quality



[primedENDO.com](http://primedENDO.com)



# MultiCROC® Advantage

**MultiCROC® Single Use Biopsy Forceps** are engineered & manufactured to collect five biopsy samples of the upper or lower digestive tract in one pass through the channel of the endoscope.

The devices are supplied sterile and are intended for single use only.

Consider the following features and benefits of the **MultiCROC®** and discover why it is the finest multi-sampling biopsy forceps available.

## FEATURES & BENEFITS

- **Five samples in one pass through the endoscope**  
Reduced procedure time benefits the patient and the provider  
Quick and accurate reach of the biopsy sites assures better diagnostics
- **Alligator style 2.4mm diameter jaws**  
Excellent tissue sampling for accurate histopathological diagnosis  
Minimum channel size of 2.8mm
- **Patented design with large biting force**  
Superior tissue quality with well-defined borders and minimal crush artifact
- **Unique geometry of the jaws**  
Easy transfer of specimens to fixative

PART NUMBER	JAWS DIAMETER	LENGTH	UNITS PER BOX
<b>Gastroscope</b>			
MC624A-20	2.4mm	160cm	20
<b>Colonoscope</b>			
MC924A-20	2.4mm	230cm	20

Read Instructions for Use Pamphlet before using.

### ORDERING INFORMATION

- ✉ Email: [sales@primedENDO.com](mailto:sales@primedENDO.com)
- ☎ Phone Toll Free: **877-565-0565**
- ✉ Fax Toll Free: **877-565-0566**

Manufactured by:  
 PriMed Instruments Inc.  
 1080 Tristar Drive, Unit #14  
 Mississauga, ON L5T 1P1  
 Canada

☎ **905-565-0565**  
 ✉ **905-565-0566**



Patents Pending / ©2019 PriMed Instruments Inc. Manufactured under ISO 13485 & MDSAP Quality Systems. CE Mark, Health Canada & FDA listed.